

Technical Bulletin

Thermolabel

Temperature Indicating Tape

Product Description:


Thermolabels are self-adhering colour changing labels which measure the surface temperature of dishware in dishwashing machines.

Using Thermolabels on a dish surface with each wash cycle will indicate the proper sanitation temperature has been reached.

Application:

- Completely waterproof
- Label changes from white to black when its calibrated temperature is reached.
- Accurate — within 1% of calibrated temperature
- Temperatures available — 71°C, 77°C, 82°C
- Large, very visibly labels with adhesive backing
- Easy to use
- Used labels can be kept for record keeping and documentation
- Low per test cost

The Thermolabel is attached to a clean dry dish and then run with a load of soiled dishes through a regular wash cycle. The label changes from white to black when the indicated temperature is reached. This is a permanent colour change, indicating the dish surface has reached the label temperature.



4 Labels per Sheet, 4 Sheets per Pack
Code: TL3160170180

Label before sanitisation

160° F	71° C
170° F	77° C
180° F	82° C

Label after sanitisation

160° F	71° C
170° F	77° C
180° F	82° C

SHELF LIFE: As a quality assured manufacturer, Castle Chemicals has a stringent Quality assurance programme. As part of this regime, the label on this product shows a batch number and date of manufacture. This product has a shelf life of 24 months from the label printed date of manufacture. This information contained herein is based on data considered accurate. However, no warranty is expressed or implied regarding the accuracy of these data or the results to be obtained from the use thereof. Castle Chemicals assumes no responsibility for personal injury or property damage to vendees, users or third parties caused by the material. Such vendees or users assume all risks associated with the use of material. **Page 1 of 1**